

1 2 3 4 5 6 7 8

REV	ECN	DESCRIPTION	DESIGN	DATE

A


A

B

B

C

C



技术要求

- 1、材质：白卡纸300g；
- 2、外表面需哑膜处理；
- 3、文字大小、颜色均按电子图档要求；
- 4、外尺寸为150*70*25mm，若图纸尺寸有问题，请供应商提供外尺寸调整；

THIS DRAWING AND THE INFORMATION CONTAINED HEREIN IS THE PROPERTY OF Singuway Biotech Inc., AND SHALL NOT BE USED, REPRODUCED, PUBLISHED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN CONSENT OF Singuway Biotech Inc.

UNIT	mm	Singuway Biotech Inc.			
APPROVALS	DATE	TITLE	单人份抗原试剂盒 (德语)		
DESIGN		DWG No	C-401-000080-01		
CHECK		MAT	白卡纸	REV	1.0
APPROVED		SCALE	1:1	SIZE	A4
SHEET 1 OF 1			PROJECTION	第一角画法	ATTR
					包材

D

1 2 3 4 5 6 7 8